

Western Blot 1 X

Western Blot 2 X

Western Blot 3 X

+

5

PERFORM
WESTERN
BLOT


1. Sample Prep 2. Prepare Gel 3. Load Gel 4. Run 5. Transfer 6. Blot 7. Develop

Samples

1. Sample Name + Lysate Type
2. Sample Name + Lysate Type
3. Sample Name + Lysate Type
4. Sample Name + Lysate Type
5. Sample Name + Lysate Type
6. Sample Name + Lysate Type
7. Sample Name + Lysate Type
8. Sample Name + Lysate Type
9. Sample Name + Lysate Type
10. Sample Name + Lysate Type
11. Sample Name + Lysate Type
12. Sample Name + Lysate Type
13. Sample Name + Lysate Type
14. Sample Name + Lysate Type
15. Protein Ladder

Protein X X

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15



250
150
100
69
50
37
25
20
15
10

BLOTTING CONDITIONS

Primary Antibody:
Goat anti-Protein X

Secondary Antibody:
Rabbit anti-Goat

Exposure Time

1 sec 1 hour

RE-PROBE



◀ SAMPLE PREP